

## 9.2 Probe List

---

The ultrasound image scanner uses probes to obtain graphic data of the human body and then displays it on the screen. Always use application-specific probes in order to obtain the best quality images. It is also important to configure the probe with the best settings for the particular organ being scanned

### 9.2.1 Probe Application and Preset

Probes, applications and settings available for this product are as follows.

Probe	Application
C2-4/20	Abdomen, Cardiac
CN4-9	Neonatal, Pediatric, Vascular
CN2-8	Abdomen, OB, Gynecology, Renal, Fetal Heart
LN5-12/40	MSK, Small Parts, Vascular
L5-12/60	Abdomen, Cardiac, Pediatric Cardiology, TCD
LV5-12/60	OB, MSK
LV2-5/180	MSK
EC4-9	OB, Gynecology, Urology

#### NOTE

---

- Besides, user can have or select own preset between User 1- 5.
- Refer to 'Chapter 3. Starting Diagnosis' for how to set probe and its preset.

### 9.3 Thermal Index (TI Table) Tables

TI(TI;Thermal Index) values displayed on the screen title bar can change depending on probes and applications. EKO7 decides automatically which TI value will be displayed out of TIs(TI;Thermal Index System), Tib(TI;Thermal Index Bone), and TIc(TI;Thermal Index Cranium). The TI values are as follows.

Probe	Application												
	Obstetrics	Gynecology	Abdomen	Fetal Heart	Cardiac	Breast	Vascular	Urology	Renal	Small Parts	Neonatal	Pediatric	Musculoskeletal
C2-4/20			TIs		TIs								
CN4-9							TIs				TIc	TIs	
CN2-8	Tib	TIs	TIs	Tib									
LN5-12/40						TIs	TIs			TIs			TIs
L5-12/60						TIs	TIs			TIs			TIs
LV5-12/60	Tib												TIs
LV2-5/180													TIs
EC4-9	Tib	TIs						TIs					